

FIG. 1

2/15




FIG. 2




FIG. 3

3/15




FIG. 4

	VT-1 ZONE	VT ZONE	VF ZONE
1-ZONE CONFIGURATION	NONE		
2-ZONE CONFIGURATION		V RATE > A RATE AFIB RATE THRESHOLD STABILITY ONSET SRD SHOCK IF UNSTABLE <sup>a</sup>	NONE
3-ZONE CONFIGURATION	V RATE > A RATE AFIB RATE THRESHOLD STABILITY ONSET SRD	SHOCK IF UNSTABLE	NONE

<sup>a</sup>. SHOCK IF UNSTABLE CANNOT BE PROGRAMMED ON IN THE SAME ZONE AS OTHER DETECTION ENHANCEMENTS THAT ARE PROGRAMMED TO INHIBIT THERAPY (ONSET, STABILITY, AND AFIB RATE THRESHOLD)

FIG. 5

5/15




FIG. 6

DETECTED RHYTHM	AFIB RATE THRESHOLD WITH STABILITY
GRADUAL, UNSTABLE A > AFIB RATE THRESHOLD	INHIBIT UNTIL V > A OR SRD (IF ON)
GRADUAL, UNSTABLE A < AFIB RATE THRESHOLD	INHIBIT UNTIL V > A OR SRD (IF ON)
SUDDEN, UNSTABLE A > AFIB RATE THRESHOLD	INHIBIT UNTIL STABLE, V > A OR SRD (IF ON)
SUDDEN, UNSTABLE A < AFIB RATE THRESHOLD	TREAT IMMEDIATELY AT END OF DURATION <sup>a</sup>
GRADUAL, STABLE A > AFIB RATE THRESHOLD	TREAT IMMEDIATELY AT END OF DURATION <sup>b</sup>
GRADUAL, STABLE A < AFIB RATE THRESHOLD	INHIBIT UNTIL V > A OR SRD (IF ON)
SUDDEN, STABLE A > AFIB RATE THRESHOLD	TREAT IMMEDIATELY AT END OF DURATION
SUDDEN, STABLE A < AFIB RATE THRESHOLD	TREAT IMMEDIATELY AT END OF DURATION

- a. IF V RATE > A RATE IS PROGRAMMED ON AND IS FALSE, THEN THERAPY WILL BE INHIBITED IF THE RHYTHM IS UNSTABLE
- b. IF V RATE > A RATE IS PROGRAMMED ON AND IS FALSE, THERAPY WILL STILL BE INITIATED BECAUSE THE RHYTHM IS DECLARED STABLE

FIG. 7

DETECTED RHYTHM	POSSIBLE COMBINATIONS	POSSIBLE COMBINATIONS
	ONSET AND STABILITY <sup>a</sup>	ONSET OR STABILITY
GRADUAL, UNSTABLE	INHIBIT UNTIL $V > A$ OR SRD (IF ON)	INHIBIT UNTIL RHYTHM BECOMES STABLE, $V > A$ , OR SRD (IF ON)
GRADUAL, STABLE	INHIBIT UNTIL $V > A$ OR SRD (IF ON)	TREAT IMMEDIATELY AT END OF DURATION
SUDDEN, UNSTABLE	INHIBIT UNTIL RHYTHM BECOMES STABLE, OR $V > A$ , OR SRD (IF ON)	TREAT IMMEDIATELY AT END OF DURATION
SUDDEN, STABLE	TREAT IMMEDIATELY AT END OF DURATION	TREAT IMMEDIATELY AT END OF DURATION

a. THE AND COMBINATION IS THE NOMINAL SETTING WHEN BOTH ARE ENABLED

FIG. 8




FIG. 9

8/15



9/15

92

**UT-1 Detection Enhancements**

Initial	Present	Change	Redefinition	Present	Change
V Rate > A Rate	On		Post-Shock	On	
A Ftb Rate Threshold	Off		V Rate > A Rate	10	
Stability	10	bpm	A Ftb Rate Threshold	9	bpm
And		ms	Stability	0:30	ms
Orset	9%	ms	Sustained Rate Duration		ms
Sustained Rate Duration	0:30				

FIG. 10

	RHYTHM DISCRIMINATION AVAILABLE PER ZONE		
	VT-1 ZONE	VT ZONE	VF ZONE
3-ZONE CONFIGURATION	ATRIAL TACHYARRHYTHMIA SINUS TACHYCARDIA	POLYMORPHIC VT	NONE
2-ZONE CONFIGURATION		ATRIAL TACHYARRHYTHMIA SINUS TACHYCARDIA OR POLYMORPHIC VT	NONE

FIG. 11

PARAMETER	ATRIAL TACHYARRHYTHMIA	SINUS TACHYCARDIA	POLYMORPHIC VT
V > A	ON	ON	
A FIB RATE THRESHOLD	200 MIN <sup>-1</sup>		
STABILITY (INHIBIT)	10 ms		
ONSET (INITIAL DETECTION ONLY)		9 %	
SRD INITIAL REDETECTION	0:30 MINUTES:SECONDS 0:15 MINUTES:SECONDS	0:30 MINUTES:SECONDS	
SHOCK IF UNSTABLE			30 MS

FIG. 12

12/15

PARAMETER	ATRIAL TACHYARRHYTHMIA	SINUS TACHYCARDIA	POLYMORPHIC VT
STABILITY (INHIBIT)	30 ms		
ONSET (INITIAL DETECTION ONLY)		9 %	
SRD INITIAL REDETECTION	0:30 MINUTES:SECONDS 0:15 MINUTES:SECONDS	0:30 MINUTES:SECONDS	
SHOCK IF UNSTABLE			30 MS

FIG. 13

13/15

92

UT-1 Detection

Initial	Present	Change	Redetection	Present	Change
Rate	160	145 bpm	Redetection Duration	1.0	sec
Interval	375	415 ms	Post-shock Duration	1.0	sec
Duration	2.5	sec			

☒ Detection Enhancements

☒ Atrial Tachycardia Discrimination

☒ Sinus Tachycardia Discrimination

Present On

Change

Cancel Changes

FIG. 14

14/15

Utilities Tachy Mode Monitor only

Surface Atrial Vent

AS VS AS VS AS VS AS VS AS VS

DDDR 60-90 ppm MTR 80 ppm

# Zones 1 2 3

125 bpm VT-1 2.5 sec SRD V>A A Fib. Onset ATP 5J 10J 31J AT PX2 5J 10J 31J 10J 20J 31J+

165 bpm VT 1.0 sec

Initial Detection Rate Interval Duration

88

Present Change

145 bpm 414 ms 2.5 sec

Redetection Duration Post-shock Duration

2.5 sec 1.0 sec

Present Change

1.0 sec 1.0 sec

94

90

☒ Detection Enhancements On

☒ Polymorphic VT Discrimination

Cancel Changes

System Summary Quick Check Tachy Parameters Brady Parameters Setup Therapy History Diagnostic Evaluation EP Test

FIG. 15

92

VT Detection Enhancements

Initial Shock If Unstable	Present 10	Change ms
------------------------------	---------------	--------------

Cancel Changes

Close

FIG. 16